

What is recurrent anterior urethral stricture disease?

A recurrent anterior urethral stricture is a narrowing of the urethra caused by scarring, that obstructs the flow of urine. If left untreated, the obstruction caused by a recurrent anterior urethral stricture can lead to a reduced quality of life, bladder and kidney damage, and urinary tract infections in men.³

Common symptoms include:

- Hesitancy when waiting for the stream of urine to begin
- Weak or irregular flow
- Straining to urinate
- Painful urination
- A sense of incomplete emptying
- Dribbling after urination

Ask your urologist if Optilume[®] may be the right solution for your recurrent anterior urethral stricture disease

Learn more at www.optilume.com

- 1 Steenkamp JW, Heyns CF, de Kock MLS. Internal Urethrotomy Versus Dilation as Treatment for Male Urethral Strictures: A Prospective, Randomized Comparison. *Journal of Urology*, Vol. 157:96-101. January 1997
- 2 Elliott SP, Coutinho K, Robertson KJ, D'Anna R, Chevli K, Carrier S, Aube-Peterkin M, Cantrill CH, Ehlert MJ, Te AE, Dann J, DeLong JM, Brandes SB, Hagedorn JC, Levin R, Schlaifer A, DeSouza E, DiMarco D, Erickson BA, Natale R, Husmann DA, Morey A, Olsson C and Virasoro R. One-Year Results for the ROBUST III Randomized Controlled Trial Evaluating the Optilume Drug-Coated Balloon for Anterior Urethral Strictures. *The Journal of Urology*[®] (2021), doi: 10.1097/JU.0000000000002346
- 3 Tritschler S, Roosen A, Füllhase C, Stief CG, Rübber H. Urethral stricture: etiology, investigation and treatments. *Dtsch Arztebl Int.* 2013;110(13):220-226. doi:10.3238/arztebl.2013.0220
- 4 Rosenfield K, Jaff MR, White CJ, et al. Trial of Paclitaxel-Coated Balloon for Femoropopliteal Artery Disease. *NEJM* 2015;372(2):145-53
- 5 Data on file: TP1041
- 6 Elliott SP, Virasoro R, Estrella R, et al. MP56-06 The Optilume Drug Coated Balloon for Recurrent Anterior Urethral Strictures: ROBUST I Clinical Study 3-year follow up, N = 43. *J Urol* 2021;206(3S):e971.



This document is for use in the European Union only.
Manufactured by Urotronic, Inc.


©2021 Laborie. All rights reserved.

Laborie⁷

USA

T +1 802 857 1300
+1 800 522 6743

E optilume@laborie.com
W optilume.com



Learn more

MKT-00554 [A] 01/22

Laborie⁷
FOR DIGNITY. FOR LIFE.



UROLOGY & UROGYNECOLOGY

Are you ready to break the recurrent anterior urethral stricture cycle?



OPTILUME

Optilume[®] urethral drug coated balloon



At last, a breakthrough treatment for strictures

Current endoscopic anterior urethral stricture treatments are common practice and often ineffective, leading to multiple repeat procedures with diminishing results.¹

Optilume[®] is an effective alternative to traditional endoscopic stricture treatments, offering minimally invasive, immediate relief.² Optilume[®] works by inhibiting new scar tissue growth that often recurs after dilations or direct visual internal urethrotomy (DVIU).


All this in a simple out-patient procedure.

What is Optilume[®] and how is it used?


- Optilume[®] is a small, cylindrical balloon coated with medication engineered to simultaneously open the urethra and slow or stop stricture regrowth
- Similar balloons have successfully prevented re-treatment in cardiovascular procedures for several years⁴
- Recovery is mild⁵
- No impact on erectile function²




Insertion



Dilation



Drug release




Durability⁴

Proven results

- Paclitaxel, the drug that coats the balloon, is a proven antiproliferative drug that has been effectively used in millions of patients worldwide since the early 1990's, in various applications

3-year freedom from reintervention rate

77%



▲ At 3 years,⁶ durability continued with 77% FFRI, with a 176% increase in Qmax (Baseline 5.1, 3-year 14.1) and a 65% decrease in IPSS (Baseline 25.2, 3-year 8.8). Robust I clinical study, n=43

Simple. Safe. Durable.

This changes everything